

# Heat shock transformation

Date: 2022-08-17

Tags: 3\_Cloning 1\_Wetlab

Created by: Stefanie Brands

1 / 1

( \_Written by Stefanie Brands\_ )

( \_Last update: 2022.08.18\_ )

{20.04.2022|Date of experiment}

A tube of {100| $\mu$ L|competent DH5 $\alpha$ |cell suspension} cells was thawed on ice for {10|min|incubation on ice}. {2| $\mu$ L|vector DNA|volume} containing {175|ng|(of:) plasmid|DNA} were added to the cell mixture. The tube was {flick(ed):|process} carefully 4-5 times to mix cells and DNA. The mixture was placed on ice for {30|min|incubation on ice}. A heat shock was applied at exactly {42 °C for 45 seconds|heat shock}. The transformed cells were placed on ice for {2|min|incubation on ice}. The cell suspension was filled up {to 1 mL|with room temperature SOC media|recovery media}. The tube was shaken vigorously at {37 °C for 45 min|recovery incubation}. {100| $\mu$ L|plated cell suspension} of each dilution and the {remaining (:pelleted:) cells|plated cell pellet} were spread onto an {LB-Kan|:selection plate:} and incubated {overnight at 37 °C|incubation}.

The transformation resulted in the formation of {more than 300|CFU|colony number}. {3|random clones|number of candidate clones} were transferred from the transformation plate to a {5 mL|:overnight preculture:} for minipreps and sequencing.